



**Operating Room Product Catalog**  
RF Surgical Systems

**CONSOLE, WAND, WAND COVER**

RF SURGICAL Catalog #	Description	Carton Quantity	Cartons/Case	Case Quantity
01-0002	RFS Detection Console	1	N/A	N/A
02-0029	RFS Blair Port Wand (Sterile)	10	2	20
01-0020	RFS Wand Cover (Sterile)	N/A	1	50



**STERILE COTTON PRODUCT LINE**


RF SURGICAL Catalog #	Description	Pack Quantity	Packs/Case	Case Quantity
<b>LAPAROTOMY SPONGES</b>				
L0418-04P01C-1	Lap Sponge, XRAY, RFD, 4x18, STRL, LF, DISP	5	40pk/cs	200
L1212-04P01C-1	Lap Sponge, XRAY, RFD, 12x12, STRL, LF, DISP	5	40pk/cs	200
L1818-04P01C-1	Lap Sponge, XRAY, RFD, 18x18, STRL, LF, DISP	5	40pk/cs	200
L1818-04P01C-R	Lap Sponge, XRAY, RFD, 18x18, Ring, STRL, LF, DISP	5	40pk/cs	200
L1818-04P01C-L	Lap Sponge, XRAY, RFD, 18x18, Noloop, STRL, LF, DISP	5	40pk/cs	200
L1818-04P02C-1	Lap Sponge, XRAY, RFD, 18x18, 10 PK, STRL, LF, DISP	10	20pk/cs	200
L0836-04P01C-1	Lap Sponge, XRAY, RFD, 8x36, STRL, LF, DISP	5	40pk/cs	200
L1836-04P01C-1	Lap Sponge, XRAY, RFD, 18x36, STRL, LF, DISP	5	20pk/cs	100
<b>SURGICAL GAUZE</b>				
G0404-16P01C-1	Gauze, XRAY, RFD, 4x4, 16-ply, STRL, LF, DISP	10	128pk/cs	1280
G0404-32P01C-1	Gauze, XRAY, RFD, 4x4, 32-ply, STRL, LF, DISP	10	128pk/cs	1280
G0804-12P02C-1	Gauze, XRAY, RFD, 8x4, 12-ply, STRL, LF, DISP	10	80pk/cs	800
G0804-24P02C-1	Gauze, XRAY, RFD, 8x4, 24-ply, STRL, LF, DISP	10	80pk/cs	800
G0804-32P02C-1	Gauze, XRAY, RFD, 8x4, 32-ply, STRL, LF, DISP	10	80pk/cs	800
<b>O.R. TOWELS</b>				
TOR-0101C-1	OR Towel, XRAY, RFD, Blue, STRL, LF, DISP	4	40pk/cs	160
TORG-0101C-1	OR Towel, XRAY, RFD, Green, STRL, LF, DISP	4	40pk/cs	160
TORW-0101C-1	OR Towel, XRAY, RFD, White, STRL, LF, DISP	4	40pk/cs	160



**NON-STERILE COTTON PRODUCT LINE (Bulk Packaging)**

RF SURGICAL Catalog #	Description	Pack Quantity	Case Quantity
<b>LAPAROTOMY SPONGES</b>			
L0418-04P01CN-1	Lap Sponge, XRAY, RFD, 4x18, NSTRL, LF, DISP	5	2240
L1212-04P01CN-1	Lap Sponge, XRAY, RFD, 12x12, NSTRL, LF, DISP	5	1680
L1818-04P01CN-1	Lap Sponge, XRAY, RFD, 18x18, NSTRL, LF, DISP	5	540
L1818-04P01CN-R	Lap Sponge, XRAY, RFD, 18x18, Ring, NSTRL, LF, DISP	5	540
L1818-04P01CN-L	Lap Sponge, XRAY, RFD, 18x18, Noloop, NSTRL, LF, DISP	5	540
L0836-04P01CN-1	Lap Sponge, XRAY, RFD, 8x36, NSTRL, LF, DISP	5	800
L1836-04P01CN-1	Lap Sponge, XRAY, RFD, 18x36, NSTRL, LF, DISP	5	400
<b>SURGICAL GAUZE</b>			
G0404-16P01CN-1	Gauze, XRAY, RFD, 4x4, 16-ply, NSTRL, LF, DISP	10	2800
G0404-32P01CN-1	Gauze, XRAY, RFD, 4x4, 32-ply, NSTRL, LF, DISP	10	1600
G0804-12P02CN-1	Gauze, XRAY, RFD, 8x4, 12-ply, NSTRL, LF, DISP	10	2800
G0804-32P02CN-1	Gauze, XRAY, RFD, 8x4, 32-ply, NSTRL, LF, DISP	10	800
<b>O.R. TOWELS</b>			
TOR-0101CN-1	OR Towel, XRAY, RFD, Blue, NSTRL, LF, DISP	no pack	400
TORG-0101CN-1	OR Towel, XRAY, RFD, Green, NSTRL, LF, DISP	no pack	400
TORW-0101CN-1	OR Towel, XRAY, RFD, White, NSTRL, LF, DISP	no pack	400

**L 1818 -04P 01 C N -1**



XRAY: X-Ray Detectable  
RFD: RF Detectable  
RING: Plastic Ring  
LF: Latex-Free  
NSTRL: Non-Sterile  
DISP: Disposable